

ISO13997
standard blade

FEATHER®

KAKEN Blade No. 99727


Blade for determination of resistance to cutting by sharp objects.

This blade is referenced on
international standard, ISO13997:2023.



100 blades
per box

Blades for testing protective clothing, gloves, etc.
Can be used for testing machines TDM-100, STM610, etc.



For inquiries:



FEATHER SAFETY RAZOR CO., LTD.
OVERSEAS TRADE DIVISION

3-70, OHYODO MINAMI 3-CHOME, KITA-KU, OSAKA 531-0075, JAPAN
PHONE: +81-6-6458-1638 FAX: +81-6-6458-1611
URL www.feather.co.jp/en/ E-mail overseas@feather.co.jp



9301 Cutting Tool
13455 Medical Products