



FEATHER BIOPSY PUNCH

MADE IN JAPAN


Feather Biopsy Punch offers excellent sharpness by ultra-precision sharpening technology which originates in surgical blade manufacturing.



Packaged with safety cover


- Sterilized by Gamma radiation
- Single use only

Anti-roll hexagon design



Completely-sealed package by special film

To open, peel the package off by pulling apart the two flaps as indicated.



Product No.	Inner Diameter		Package
FBP-15	1.5mm		20pcs per Box
FBP-30	3mm		
FBP-40	4mm		
FBP-50	5mm		
FBP-60	6mm		
Actual size			



FEATHER SAFETY RAZOR CO., LTD.
OVERSEAS TRADE DIVISION

3-70, OHYODO MINAMI 3-CHOME, KITA-KU, OSAKA 531-0075, JAPAN
PHONE: +81-6-6458-1638 FAX: +81-6-6458-1611
URL www.feather.co.jp/en/ E-mail overseas@feather.co.jp



3001-Cutting Tool
13485-Medical Products