

List of mountable microtomes for blade holders

In some cases FEATHER blade holders can not be fixed in older models of knife holders due to the thickness and the space at the knife holder clamping systems.

Type	Manufacturer	Model	Microtome Holder Set No.								
			160/160E	240/S/E	240/D/E	130A/AE	160C	130C	130C-S	130C-AO	
Sliding Microtome	ANGLIA	AS 200	○	○	○						
	A.O	860	○	○	○						
	ERMA	GE 68	○	○	○						
	LEITZ	1208 / 1400 STANDARD SLEDGE (Leica SM 2400)	○	○	○						
	PFM AG	Slide 2002	○								
		Slide 2003	*1 ○	*1 ○	*1 ○						
		Slide 4003E	○	○	○						
	SAKURA	IVS 410 / CRM 440	○	○	○						
		REM 700	○	○	○						
	YAMATO	ROM 380	○	○	○						
		ROM 670e	○	○	○						
		SM 2000R / SM 2010 R	○	*2 ○	○						
	LEICA / REICHERT JUNG	SM 2400 / 2500	○	○	○						
		HN 40	○								
		HM 430 / HM 450	○	○	○						
	MICROM / THERMO	HM 400S / HM 400R	○	○	○						
		HM 440E	○	○	○						
		Hacker-Bright 8000 Sledge	○	○	○						
	BRIGHT	HYRAX 5 30 / S 50	○	○	○						
	ZEISS	HYRAX AS 200 / AS 300 / AS 325 / AS 400 / AS 500	○								
Rotary Microtome	ANGLIA	AS 200 / AS 300 / AS 325 / AS 400 / AS 500	○								
	A.O	820	○			○					
	ERMA	LARGE / SMALL	*2 ○								
	LEITZ	1512 / 1516 / 1715	*2 ○								
	LIPSHAW	45	○								
	LKB	2218 / 2230	*2 ○								
	MICRO-TEC	CUT 4055 / 4050 / 4060	○								
	PFM AG	Rotary 3003	○								
		Rotary 3002	○								
	RMC	MT 910 / 920 / 940	○								
	SAKURA	Accu-Cut SRM 200	○								
		Accu-Cut SRM 400	○								
	SLEE	B / BR1 / BR R1 / BR M	○								
		Cut 4062 / 5062 / 6062	○								
	YAMATO	CYLINDER / R 240	○								
		800 / 820 / 2025 / 2035	○								
	LEICA / REICHERT JUNG	1140 / 1150 / 1165	○								
		2030 / 2040 / 2045 / 2050 / 2055 / 2065	○								
		RM 2125 / 2135 / 2145 / 2165	○								
	MICROM / THERMO	HM 315 / 320 / 325 / 330 / 340 / 350 / 360	○								
Erogostar		○									
Shandon Finesse 325 / Shandon Finesse ME+ Microtome		○									
MEDITE	Hypercut	○									
MEDITE	Meditome	○									
DIAPATH	Semiautomatic CUT R / CUT X	○									
ZEISS	HYRAX M 25 / M40 / M50 / M60	○									
Cryostat Microtome	ANGLIA	AS 600 / AS 650							○		
	A.O	CRYOCUT 2, HISTOSTAT							○	○	
	BRIGHT	PEARSE-TYPE (FS / FCS / CR)							*1 ○		
		AO BUILT-IN TYPE (FS / FAS)									○
		5030									
	DEMON / IEC	MINOTOME							○		
	LEITZ	1720 / 1720 DIGITAL							○		
	LIPSHAW	CRYOTOME							○		
	MILES	LAB TEC 2							○		
	SAKURA	COLDTOME								○	
		CM / MC 2								○	
		Tissue-Tek Cryo3								○	
	SLEE	MSK / MTE / MHR / MNT / MTC / MEU / MEV								○	
		MNT / MEV								○	
	YAMATO	CR 502								○	
		CM 1100 / CM 1510 S / CM 1850								○	
	LEICA / REICHERT JUNG	CM 1850 UV / CM 1950 / CM 3050 S								○	
		Frigocut 1800 / 2700 / 2800 / 2800E								○	
		Cryocut 3000								○	
	MICROM / THERMO / SHANDON	HM 500 / HM 525 / HM 550 / HM 560								○	
620 / 620E / 620ME / 620 SME									○		
AS 600 / AS 650									○		
MEDITE	Medistat							○			
DIAPATH	Cryostat Cryopath							○			
HISTO LINE	MC 4000							○			
ZEISS	HYRAX C25 / C50 / C60							○			
LKB	Cryostat							○			
Freezing Microtome (Sartorius Type)		FOR LARGE-SIZED FROZEN TISSUES	○		○			○		○	
		FOR SMALL-SIZED FROZEN TISSUES			○			○		○	

*1. In some cases, FEATHER blade holders cannot be mounted on older models of microtomes due to the thickness of blade fixing plate.

*2. No.160E and 240E Microtome Holder Sets are not applicable to rotary type microtomes.