



Surgical blade and handle

Blade remover

Surgery

Blade pliers

Disposable scalpel

Safeshield Scalpel

Blood lancet

Disposable dermatome

Skin graft knife blade

Biopsy punch

Skin scraper

Chokkin

Disposable prep razor

MADE IN JAPAN

Chokkin (with strap)

Pocket-sized scissors designed for exceptional handling ease, cutting ability and safety. "Chokkin" scissors are small enough to fit in a nurse's pocket. And come with a holder and clip for easy access without damage to the pocket. The perfect tool for the day-to-day hospital routine of nurses. Blade edges are specially processed with fluoroplastic. Residual glue from adhesive tape will not cling to blade after cutting, thus assuring smooth cutting reliability of the scissors.

Improved Safety


The blade tip is rounded and protected with a plastic guard. When cutting bandages, only the plastic guard touches the skin. Thus, adding relief and reassurance to the patient during cutting procedures.

Friendly Design

"Chokkin" are designed with gentle, curved lines.



Elastic strap provided



*Stainless steel blade with plastic handles (ABS resin, heatproof temperature : 70°C).

Disposable prep razor 121 S-B

*Quick and easy to use. *Individually blister packed. *Superbly sharp double-edge blade.

		Disposable prep razor 121S-B
Material	Blade	Stainless steel
	Handle	Plastic
Packaging		Supplied in box of 50 pcs.



FEATHER SAFETY RAZOR CO., LTD.
OVERSEAS TRADE DIVISION

3-70, OHYODO MINAMI 3-CHOME, KITA-KU, OSAKA 531-0075, JAPAN
PHONE:+81-6-6458-1638 FAX:+81-6-6458-1611
URL www.feather.co.jp/en/ E-mail overseas@feather.co.jp



9001 Cutting Tool
13485 Medical Products

